
## Preclinical Oncolysis Prototype

## Overview

The Openwater Preclinical Oncolysis
Prototype is designed to help researchers
investigate the effect of different ultrasound
parameters on a variety of *in vitro* and
preclinical *in vivo* targets. Certain acoustic
parameters may be well suited for
damaging cancer cells while sparing
surrounding healthy tissue, and this system
is designed to systematically explore such
effects in a laboratory setting.

Multiple focused ultrasound transducers are connected to a programmable signal generator through an amplifier and a switch to select which transducer is active. Samples are placed beneath the transducer(s), using 3D-printed coupling cones filled with water to make a path for the ultrasound and set the depth of the focal spot within the target. This system is compatible with in vitro targets held in standard 6-well plates as well as with mouse flanks, through the assistance of a stereotactic frame.



Parameter	Min	Мах	Units
Frequency	70	1000	kHz
Output Voltage	0	843	V
Peak Negative Pressure	0	1*	MPa
Mechanical Index	0	3.8*	
ISPPA	0	100*	W/cm2
ISPTA	0	1000*	mW/cm2
Burst Length	2	40	ms
Duty Cycle	0.5	10	%
Treatment Duration	5	60 * 15	seconds

\*These are the maximum values that can be set when using the parameter as the target. If the parameter is derived and not used as the target, it may exceed this range to reach another parameter's target value.

Controller software allows for configuration of

treatment parameters through a simple GUI application. With its default settings, the system can be configured to output focused ultrasound from the parameters shown in the table to the right. This allows for investigation of treatment parameters both within and well above the safety levels used for diagnostic ultrasound imaging.

## **Getting Started**

Additional information can be found in the following repositories:

- opw oncolysis hw assembly/operation instructions, mechanical drawings, 3D models
- opw oncolysis sw software to configure and operate the system
- opw oncolysis data research data and experimental procedures

CAUTION - Investigational device. Limited by Federal (or United States) law to investigational use. This system has not been evaluated by the FDA and is not designed for the treatment or diagnosis of any disease. It is provided AS-IS, with no warranties. User assumes all liability and responsibility for identifying and mitigating risks associated with use